

# **DRAGEN TSO500 Analysis Software Release Notes**

## **V1.1.1**

***For TruSight Oncology 500 Assay***

**April 13, 2022**

## Introduction

These Release Notes detail the key features and known limitations to software components for the DRAGEN TSO500 v1.1.1 Analysis Software.

This software is intended for use with the TruSight Oncology 500 Assay.

- Software Version: 1.1.1
- Docker Image ID: f3e439815738

The software installer script, *install\_DRAGEN\_TSO500-1.1.1.run*, includes the following:

- *dragen\_tso500\_1.1.1.tar* – a tar file of the DRAGEN TSO500 docker image.
- *uninstall\_DRAGEN\_TSO500-1.1.1.sh* – a script for uninstalling DRAGEN TSO500.
- *check\_DRAGEN\_TSO500-1.1.1.sh* – a script for self-testing DRAGEN TSO500.
- *build-hashtable\_DRAGEN\_TSO500-1.1.1.sh* – a script for building the hash table.
- *install.sh* – a script used to install DRAGEN TSO500 based on the contents listed.
- *resources/* – a directory containing all resources files necessary for DRAGEN TSO500.
- *dragen-3.6.6-2.el7.x86\_64.run* – the DRAGEN installer.

### NEW FEATURES:

- Initial release

### DEFECT REPAIRS (COMPARE TO DRAGEN TSO500 ANALYSIS SOFTWARE V1.0.0 (EARLY ACCESS LIMITED RELEASE VERSION)):

- The fusion caller in the Early Access Limited Release version of the software may not always call fusions with expression of intronic region or fused in antisense orientation. The DRAGEN TSO500 v1.1.1 analysis software does not have this limitation.
- The fusion caller in the Early Access Limited Release version of the software may classify some RNA read-through transcripts with small insert size as supporting reads for a fusion in the opposite orientation if the neighboring genes are in close proximity and have local repeat sequences to seed the initial fusion candidate. The DRAGEN TSO500 v1.1.1 analysis software does not have this limitation.
- The Early Access Limited Release version of the software on Illumina Connected Analytics required output for all samples scattered to be gathered at once. Gathering a subset of results was not supported. The DRAGEN TSO500 v1.1.1 analysis software does not have this limitation.

### KNOWN ISSUES:

- Moving or modifying files during the analysis may cause the analysis to fail or provide incorrect results.

- Using control-c during a running analysis may cause an FPGA error. To recover from an FPGA error, shut down and restart the server.
- The sample sheet should not have blank rows between samples in the [Data] section. Blank rows in the [Data] section may cause a run failure.
- Performance not verified using reads other than 2 x 101.
- The TSO500 RNA workflow is unstranded. Fusions or splice variants could involve antisense transcripts instead of the reported genes.
- Fusion caller may not always call fusions where breakpoint(s) are located in region(s) with high homology.
- The cloud workflow will fail if blank rows are present after the [Data] section.

**PRODUCT LIMITATIONS:**

- The sample sheet must be configured as described in the User Guide.
- The values in the Run Metrics section will be listed as 'NA' if the analysis was started from FASTQs or if the analysis was started from BCLs but the InterOp files are missing or corrupted.
- Unmapped long insertions are not likely to occur on shorter indels because there is sufficient reference-matching sequence in the reads. Product claims only indels up to 25 base pairs.
- Complex variants are specifically output only for a specific region of the EGFR gene, component and phased variants would both be contained in the output.
- Incorrect calculation of variant allele frequency can occur in variants near the start and end of genomic reads, but there is a low probability of incorrect variant allele frequency in called variants due to sufficient variation in read start and end positions.
- Germline estimation uses latest publicly available population data and is estimated to be representative of targeted population. The impact of rare germline mutations is expected to be limited.
- Germline estimation is difficult when tumor purity is >85% causing expected variant allele frequency for somatic and germline variants to converge.
- Poor quality wild type reads may align as chimeric and be miscalled during RNA analysis.
- The Illumina Annotation Engine (aka Nirvana) may report incorrect protein (P-Dot) and transcript (C-Dot) changes in HGVS nomenclature for small variants located on a RefSeq transcript where an RNA-edit has occurred. Most known variants on these transcripts are unaffected. A list of affected Canonical RefSeq transcripts and Cosmic Variants from those transcripts is provided below. A full explanation of this product limitation is provided in PQN2020-1090. [1]

## Affected Canonical RefSeq Transcripts

<b>Transcript ID</b>	<b>Gene Symbol</b>
<b>NM_002467.4</b>	<b>MYC</b>
<b>NM_003224.5</b>	<b>ARFRP1</b>
<b>NM_004119.2</b>	<b>FLT3</b>
<b>NM_006904.6</b>	<b>PRKDC</b>
<b>NM_198291.2</b>	<b>SRC</b>
<b>NM_021960.4</b>	<b>MCL1</b>
<b>NM_001025366.2</b>	<b>VEGFA</b>

Affected Cosmic Variants from Canonical RefSeq Transcripts

The list of affected variants is based on an analysis of COSMIC database version 92 variants located along the Canonical RefSeq Transcripts listed above [2]. New variants are regularly submitted to COSMIC, and this list of affected variants may change over time.

Chr:Position	REF*	ALT**	Gene Symbol	Transcript ID	COSMIC_ID
chr1:150548890	A	ATCTA	MCL1	NM_021960.4	COSV57189597
chr6:43738444	C	T	VEGFA	NM_001025366.2	COSV104569261
chr8:48805817	G	GG	PRKD C	NM_006904.6	COSV58041377
chr8:128748839	GC	G	MYC	NM_002467.4	COSV104388447
chr8:128748840	C	A	MYC	NM_002467.4	COSV104388806
chr8:128748840	C	G	MYC	NM_002467.4	COSV104388204
chr8:128748841	T	C	MYC	NM_002467.4	COSV104388663
chr13:28608094	C	CACTTTTCCAAAAGCACCTG ATCCTAGTACCTTCCAAAC TCTAAATTTTCTCTTGAA CTCCATTTGAGATCATATT CATATTCGTTTCATC	FLT3	NM_004119.2	COSV54069050
chr13:28608124	C	CTTCCCAAACCTACTGTTG CGTTCATCACTTTTCCAAA GCACCTGATCCTAGTACC	FLT3	NM_004119.2	COSV54044227
chr13:28608129	C	CAAACCTCAAAGCACCTGAT CCTAGTACCTTCCC	FLT3	NM_004119.2	COSV54054381
chr13:28608129	C	CAAACCTCTAAATTTTCTTTG GAAACTCCCATTATCCTAGT ACCTTCCC	FLT3	NM_004119.2	COSV54043729
chr13:28608129	C	CAAACCTCTAAATTTTCTTTG GAAACTCCCATTTTCCAAA GCACCTGATCCTAGTACCTT CCC	FLT3	NM_004119.2	COSV54075746
chr20:36030939	G	GTGGCC	SRC	NM_198291.2	COSV99050886
chr20:62331336	C	CC	ARFR P1	NM_003224.5	COSV53926174

\*Reference base(s)

\*\*Alternate base(s)

[1] DRAGEN TSO500 uses the Canonical RefSeq transcript when annotating variants passed into the Combined Variant Output file. The Illumina Annotation Engine selects canonical transcripts based on the following rules:

- Order all overlapping transcripts by coding sequence length.
- Pick the longest transcript that has an associated Locus Reference Genome (LRG) sequence.
- If no LRGs exist for the set of transcripts, pick the longest transcript that is coding.
- If there is a tie, pick the transcript with the smaller accession id number.

[2] Released 27 August 2020.

## Release History

Version	ER#	Author	Description of Change
00	CN 1067391	Mario Duff	Initial Release