VeriSeq NIPT Solution v2



Consumables and Equipment List

Consumables

Consumable	Supplier
1000 µl Conductive Non-Sterile Filter Tips	Hamilton, part # 235905
300 µl Conductive Non-Sterile Filter Tips	Hamilton, part # 235903
50 µl Conductive Non-Sterile Filter Tips	Hamilton, part # 235948
Deep-well reservoir with the following specifications: • SLAS 1–2004 microplate format with 96 pyramidal or	General lab supplier
 conical bottom wells and a 240 ml minimum capacity. Polypropylene with preference for low DNA binding for all sample contact surfaces. Internal dimensions (liquid level) are compatible with automated aspiration and dispensing steps of VeriSeq NIPT Microlab STAR. Height dimensions are compatible with automated movements of VeriSeq NIPT Microlab STAR. 	Compatible reservoirs: Corning Axygen, product # RES-SW96-HP-SI Agilent, product # 201246-100
Reagent tub with the following specifications: • Tub that fits securely, but not forced, into the carrier of	General lab supplier
the Alleria Committee of the Atlene to the CTAD with the committee of the Atlene to the	Campatible tube.

- Tub that fits securely, but not forced, into the carrier of the VeriSeq NIPT Microlab STAR with tapered bottom and a 20 ml minimum capacity.
- Polypropylene that is free from RNase/DNase.
- Internal reservoir dimensions (liquid level) generate liquid levels using assay reagent volumes that are compatible with automated aspirating and dispensing steps of VeriSeq NIPT Microlab STAR.
- Height dimensions are compatible with automated movements of VeriSeq NIPT Microlab STAR.

Compatible tubs:

 Illumina Reagent Tub, product # 020095418



Consumable	Supplier
Deep-well plates with the following specifications: • SLAS 1–2004, 3–2004, and 4–2004 microplate format	General lab supplier
 SLAS 1-2004, 3-2004, and 4-2004 inicropiate format with 96 pyramidal or conical bottom wells and a 2 ml minimum well capacity. Translucent polypropylene, with preference for low DNA binding material for all sample contact surfaces. Well dimensions generate a liquid level that is compatible with automated aspiration and dispensing steps of VeriSeq NIPT Microlab STAR. Plate skirt that permits placement of plate barcodes to require position with secure, flat surface adhesion. Torque-resistant frame able to sustain a minimum of 5600 × g. Plate height dimensions are compatible with automated movements of VeriSeq NIPT Microlab STAR 	Compatible plates: • Eppendorf, part # 0030505301 • Eppendorf, part # 30502302 • USA Scientific, part # 1896-2000
384-well plate with the following specifications:	General lab supplier

- Microplate with 384 wells, optimized for low-volumes, with a 50 µl minimum well capacity.
- Black opaque polystyrene with light-blocking and low DNA binding for all sample contact surfaces.
- Well dimensions generate liquid levels that are compatible with automated aspiration and dispensing steps of VeriSeq NIPT Microlab STAR.
- Plate height dimensions are compatible with automated movements of VeriSeq NIPT Microlab STAR.
- Plate skirt that permits placement of plate barcodes to required position with secure, flat surface adhesion.

Compatible plates:

Corning, product # 3820

Consumable	Supplier
 96-well plate with the following specifications: Microplate with a torque-resistant frame able to sustain a minimum of 5600 × g and 96 translucent wells with tapered bottoms, raised rims, and a 150 µl minimum well capacity. Polypropylene that is free from RNase/DNase with low DNA binding for all sample contact surfaces. Well dimensions generate liquid levels that are compatible with automated aspirating and dispensing steps of VeriSeq NIPT Microlab STAR. Plate height dimensions are compatible with automated movements of VeriSeq NIPT Microlab STAR. Plate skirt that permits placement of plate barcodes to required position with secure, flat surface adhesion. Compatible with thermal cyclers for denaturing. 	General lab supplier Compatible plates: Eppendorf, part # 0030129512 Eppendorf, part # 30129580 Eppendorf, part # 30129598 Eppendorf, part # 30129660 Eppendorf, part # 30129679 Bio-Rad, part # HSP9601
One of the following seals: Microseal 'F' Foil Foil seals	Bio-Rad, catalog # MSF1001 Beckman Coulter, item # 538619
DNase/RNase-free water	General lab supplier
Ethanol, 100% (200 proof), molecular biology grade*	General lab supplier
Sequencing reagents and consumables required for the next-generation sequencing (NGS) system.	
If using a NextSeq 550Dx Sequencing System: NextSeq 550Dx High Output Reagent Kit v2.5, 75 cycles	Illumina, part # 20028870
Cell-Free DNA BCT CE	Streck, catalog # 218997
Push Caps	Sarstedt, order # 65.802
2 ml Screw-cap tubes	General lab supplier
20 µl filter tips for 20 µl pipettor	General lab supplier
200 µl filter tips for 200 µl pipettor	General lab supplier
1000 µl filter tips for 1000 µl pipettor	General lab supplier
Equivalent: • An alcohol-based rapid disinfectant spray • A solution of disinfecting detergent Recommended: • Deionised water and 70% ethanol	General lab supplier

^{*} Nonmolecular biology grade ethanol can negatively impact performance of the assay.

Optional Consumables

Consumable	Supplier
Dulbecco's Phosphate-Buffered Saline (DPBS) for no template control (NTC)	General lab supplier
Tube, screw cap, 10 ml (for control samples only)	Sarstedt, order # 60.551
Tube, screw cap, 50 ml	General lab supplier
25 ml serological pipettes	General lab supplier
10 ml serological pipettes	General lab supplier

Equipment Required, Not Provided

Equipment	Supplier
 A next-generation sequencing (NGS) system with the following capabilities: 2 x 36 bp paired-end sequencing Compatible with VeriSeq NIPT Sample Prep dual index adapters Automatic production of BCL files Two channel chemistry 400 million paired-end reads per run Compatible with VeriSeq NIPT Assay Software v2 or a NextSeq 550Dx Sequencing System. 	Instrument supplier or Illumina, part # 20005715
Basic lab equipment: safety glasses, lab coats, bucket	powder-free protective gloves, stopwatch or timer, ice
Freezer, -25°C to -15°C	General lab supplier
Microcentrifuge	General lab supplier
Pipette aid	General lab supplier
Refrigerator, 2°C to 8°C	General lab supplier
20 µl single-channel pipettes	General lab supplier
200 µl single-channel pipettes	General lab supplier
1000 µl single-channel pipettes	General lab supplier
Vortexer	General lab supplier



Equipment	Supplier
Centrifuge and rotor assembly for blood collect	ion tubes
 Equivalent: Refrigerated centrifuge capable of 1600 x g with no-brake option Swinging bucket rotor with buckets Bucket inserts, 24, 48, or 96 tube capacity, 76 mm minimum depth Insert adapters to support 16 x 100 mm blood collection tubes Recommended: 	General lab supplier
 Allegra X12R Series Centrifuge, 1600 g Allegra Centrifuge GH-3.8 Rotor with 	Beckman Coulter, item #392304 (120 V or 230 V)
buckets • Allegra Centrifuge Bucket Covers, set of	Beckman Coulter, item #369704
two	Beckman Coulter, item #392805
 Allegra Centrifuge Adapter Assembly, 16 mm, set of four 	Beckman Coulter, item #359150
Centrifuge and rotor assembly for microplates	
 Equivalent: Centrifuge capable of 5600 × g Swinging plate rotor with 96-well plate carriers, 76.5 mm minimum depth. 	General lab supplier
 Multifuge X4 Pro-MD 120V TX-1000BT Sorvall Legend XTR Centrifuge 	Thermo Scientific VWR catalog number 76326-256 Thermo Fisher Scientific, catalog # 75004521 (120 V) or catalog # 75004520 (230 V)
 HIGHPlate 6000 Microplate Rotor Rotor high plate 6000 Support base for microplates Recommended: 	Thermo Fisher Scientific, catalog # 75003606 Thermo Scientific VWR catalog number 97040-244
MicroAmp 96-Well Support Base96-Well PCR Plate Carrier	Thermo Fisher Scientific, catalog # 4379590 Thermo Fisher Scientific, catalog # AB-0563/1000



Equipment	Supplier
One of the following microplate readers, or equivalent, (fluorometer) with SoftMax Pro v6.2.2 or higher: Gemini XPS SpectraMax M2, M3, M4, and M5. Purple insert is included with microplate reader for use in workflow.	Molecular Devices, part # XPS Molecular Devices, part # M2, M3, M4, and M5
SpectraMax High-Speed USB, Serial Adapter	Molecular Devices, part # 9000-0938
Thermal cycler with the following specifications: • Heated lid • 4°C to 98°C temperature range • ±2°C temperature accuracy • 2°C per second minimum ramp rate • Compatible with Twin.tec PCR Plate 96-well, full skirt	General lab supplier
VeriSeq NIPT Microlab STAR	Hamilton, part # 95475-01 (115 V), part # 95475-02 (230 V), or part #806288 (for Hamilton Company Bonaduz)
VeriSeq Onsite Server v2 or an upgraded VeriSeq Onsite Server	Illumina, part # 20028403 or 20047000 (v2) or #15076164 or # 20016240 (upgraded)

Optional Equipment

Equipment	Supplier
Pluggo Decapper System	LGP Consulting, part # 4600 4450
SpectraMax SpectraTest FL1 fluorescence validation plate	Molecular Devices, part # 0200-5060
Tube Revolver/Rotator, 15 ml tubes, 40 rpm, 100–240 V	Thermo Scientific, catalog # 88881001 (US) or catalog # 88881002 (EU)

VeriSeq NIPT Sample Prep

Consumable	Part #
VeriSeq NIPT Sample Prep (24 Samples)	20025895

Consumables and Equipment List



Consumable	Part #
VeriSeq NIPT Sample Prep (48 Samples)	15066801
VeriSeq NIPT Sample Prep (96 Samples)	15066802