

NeoPrep Software Release Notes

NeoPrep System Suite v1.1.2.10

NeoPrep Protocols v1.1.2.12

NeoPrep Control Software 1.1.0.8

BaseSpace Broker 2.7.0.30

For NeoPrep Systems

Introduction

These release notes describe key changes to software components for the NeoPrep Library Prep System. If the system is already running System Suite v1.1.1.4, the System Suite v1.1.2.10 update only adds quantification and normalization (quant and norm) to the mRNA protocol and is not required if that assay will not be run, or will be run without quant and norm.

The v1.1.2.10 software package includes:

- NeoPrep System Suite 1.1.2.10
- NeoPrep Protocols v1.1.2.12
- NeoPrep Control Software 1.1.0.8
- BaseSpace Broker 2.7.0.30

Protocols Software v1.1.2.12

New Features in v1.1.2.12 (if already running v1.1.1.4)

- Adds TruSeq Stranded mRNA quant and norm functionality.