

# **Correlation Engine**

## **Customer Release Notes**

**Version 2.4.0**

**August 6, 2024**

## Introduction

These Release Notes detail the latest release of Correlation Engine, including known issues.

Correlation Engine provides an interface with tools and content for exploratory, preclinical, and early translational research in drug discovery and biotechnology research. The software is a multidimensional integration platform designed to search and correlate global collections of heterogeneous, large-scale experimental data across diverse platforms, organisms, and therapeutic areas.

### NEW FEATURES:

- Heatmap displays of gene and biogroup Meta-Analysis results
- Updates to "What's New at Correlation Engine?"

### FIXED ISSUES:

- Fixed an issue where study counts were mismatched in Disease Atlas
- Fixed an issue where some import scoring jobs failed and required re-submission
- Fixed an issue where deleted biosets still affected overall study ranks
- Fixed an where some studies with the Alzheimer's disease tag were not categorized

### KNOWN ISSUES:

- When searching by gene query on the cell lines tab in Body Atlas, some genes return no results in the expanded third level details

### DATA CONTENT UPDATE:

	Added since 6-Feb-24 release	Total public content
Studies	210	26,117
Biosets	1,094	188,977